



Invisible Products
(IIC)



Receiver In Canal
(RIC)




Completely In Canal
(CIC)



Behind The Ear
(BTE)



In The Canal
(ITC)



In The Ear
(ITE)

Style	Invisible Products (IIC)	Receiver In Canal (RIC)	Completely In Canal (CIC)	Behind The Ear (BTE)	In The Canal (ITC)	In The Ear (ITE)
Available Colours	Variety	Variety	Variety	Variety	Variety	Variety
Wireless Connectivity	✓	✓	✓	✓	✓	✓
Rechargeable Solutions		✓				
Tinnitus Solutions	✓	✓	✓	✓	✓	✓
Single-Sided Hearing		✓		✓		
Technologies Available	SoundLens Synergy iQ® SoundLens Synergy®	Livio® Livio AI® Halo iQ™ Muse iQ™ Halo 2™ Muse™ Xino Wireless™	Muse iQ™ Muse™	Livio® Livio AI® Muse iQ™ Halo™ Muse™	Muse iQ™ Muse™	Muse iQ™ Muse™
Hearing Loss	Mild to Moderately Severe	Mild to Severe	Mild to Moderate	Moderate to Severe	Mild to Mildly Severe	Mild to Severe
Battery Size	10	312 or 13	10 or 312	312 or 13	312	312 or 13